

Aspire Bath Transfer Bench - Bariatric



- Reinforced design for the bariatric user with a maximum weight of 227 kg
- It has a textured moulded plastic seat and backrest
- Fully height adjustable with a side arm for secure transfers
- Adjustable backrest to suit right or left facing tubs
- Suction cup foot pieces secure the bench to the tub floor
- Lightweight aluminium frame construction

Warranty Information

- 1 Year

✓ Contracted DVA product

✓ Contract approved

ndis

Item Code:

05_123121160_0103_1_2

Shower Support –

Composite:05_415_0103_1_2

Click [here](#) for more information



Compliance

AS/NZS 3973:2009

Code	Display Name	RRP	Box Quantity	Seat Depth	Seat Width	Seat To Floor Height	Safe Working Load	Unit Weight
BTS105300	Aspire Bath Transfer Bench - Bariatric	329	1	405 mm	670 mm	485-595 mm	227 kg	7.2 kg