

## Invacare Birdie Compact Patient Lifter



- Birdie “Compact” has a 55 cm width and can be folded and unfolded without the need for tools or special skills required. This makes it ideal for the nursing or home care environments where space and time are at a premium
- The Birdie Compact is the only lifter that is TÜV certified that is this narrow – unique to the market
- All attachments on the Birdie Compact are quick release, making it very easy to change when required
- Intelligent control system includes a service reminder, emergency operation, audio battery alarm, and battery indicator
- Manual leg width adjustment
- Safe patient handling with up to 150 kg
- Supplied with yoke attachment
- Electric emergency lowering up and down
- Item Weight: 36 kg

### Optional Accessories:

- Pivot Frame (LSS364112)
- Full range of slings

## Warranty Information

- 3 Years Electronics
- 2 Years Battery
- 1 Year Handset

- ✓ Contracted DVA product
- ✓ Contract approved

### ndis

#### Item Code:

05\_123603161\_0105\_1\_2

#### Hoist –

**Composite:**05\_420\_0105\_1\_2

Click [here](#) for more information

### Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Total Length	Total Width	Total Height	Lift Range	Leg Gait Range	Safe Working Load	Material
LSS364100	Invacare Birdie Compact Patient Lifter	1000 mm	550 mm	1600 mm	530 - 1600 mm	550 - 890 mm	150 kg	Steel