

## MOTOmed Gracile 12



**RECK**  
MOTOmed®

Effective movement therapy for children with special needs, it is the only motorised movement therapy device that fulfills the special needs of children.



- The Gracille12 comes with both upper limb and lower limb trainers installed
- During active training, children can play MOTOmax for more fun during training.
- Leg trainer (optional: arm/upper body trainer).
- Paediatric safety foot shells with paediatric leg guides.
- Small distance of 12 cm between the inner rims of the foot shells.
- Two-stage pedal radius adjustment (3.5 cm / 7.0 cm).
- Height adjustable pedal axle.
- Handlebar adjustable in height and inclination.
- Mobile with big castors (13cm diameter).
- All-metal construction, premium and stable.
- Additional stability by telescopic front stand.
- MOTOmed SmoothDriveSystem with a gentle startup and rundown.
- MovementProtector & SpasmControl.
- Speed regulation from 0 to 60 rpm (passive).
- Regulation of motor power from 1 to 10.
- Preset the therapy duration from 0 to 120 minutes.
- User-friendly operating panel with a colour coded system and big, easy-feel buttons.
- Language selection between 27 different languages.

| Code      | Display Name                                     | Total Length | Total Width | Total Height  | Safe Working Load | Unit Weight |
|-----------|--|--------------|-------------|---------------|-------------------|-------------|
| FPM101112 | MOTomed Gracile 12 - Custom (as per script)      | 60 - 88 cm   | 55 cm       | 99 - 114 cm   | 95 kg             | 33 kg       |
| FPM101110 | MOTomed Gracile w/ Arms - Custom (as per script) | 600 - 880 mm | 550 mm      | 990 - 1140 mm | 90 Kg             | 33 Kg       |