

Conni Reusable Bed Pad with Tuck-ins



- The Conni absorbent and waterproof bed pads with tuck-ins assist in the management of incontinence, bedwetting and toilet training
- Designed to hold 2.5 litres over an eight hour period, which means they are the perfect addition to any single or king-single beds that require full void protection
- Cost-effective and environmentally sustainable alternative to disposable.
- Certified as being free from harmful materials and are safe even against newborn skin
- Soft and breathable
- Washable and reusable



ndis

Item Code:

03_092489060_0103_1_1

Click [here](#) for more information

Compliance

No applicable compliance standard

Warranty Information

- 1 Year



Code	Display Name	RRP	Box Quantity	Total Length	Total Width	Colour	Capacity
CON228510	Conni Bed Pad with Tuck-ins - Mauve	79	11	1000 mm	1000 mm	Mauve	2.5 litres
CON228515	Conni Bed Pad with Tuck-ins - Teal Blue	79	11	1000 mm	1000 mm	Teal Blue	2.5 litres
CON228516	Conni Bed Pad with Tuck-ins - White	-	-	1000 mm	1000 mm	White	2.5 litres
CON228527	Conni Bed Pad with Tuck-ins - Animal print	-	-	1000 mm	1000 mm	Animal print	2.5 litres