

Aspire Deluxe General Purpose Sling - Mesh



- Australian standards tested and registered with the TGA
- Versatile sling offering body contoured fabric shaping for maximum patient comfort and support
- Quality - manufactured in the UK using European materials, computer guided stitching and CNC fabric cutting
- Durable - loop webbing has a breaking strain in excess of one tonne and is designed so that it will not tangle or “cork-screw” during cleaning
- Serial Numbered - slings are embroidered with a unique serial number to maintain traceability and compliance at all times
- Service & Maintenance - service record labels ensure you can quickly and easily identify service history
- Easy Sling Identification - all slings are colour coded by size and have reinforced edging for longevity
- Easy Clean - machine washable up to 85°C temperature
- For FIM score 1 or 2

- ✔ Contracted DVA product
- ✔ Contract approved

ndis

Item Code:

05_123603164_0105_1_2

Sling – Composite:

05_421_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Warranty Information

- 1 Year

| Code | Display Name | RRP | Box Quantity | Safe Working Load | Material |
|-----------|--|-----|--------------|-------------------|----------|
| LSS390716 | Deluxe General Purpose Sling - Mesh - Small | 355 | 1 | 200 kg | Mesh |
| LSS390717 | Deluxe General Purpose Sling - Mesh - Medium | 355 | 1 | 200 kg | Mesh |
| LSS390718 | Deluxe General Purpose Sling - Mesh - Large | 355 | 1 | 300 kg | Mesh |

| Code | Display Name | RRP | Box Quantity | Safe Working Load | Material |
|-----------|---|-----|--------------|-------------------|----------|
| LSS390719 | Deluxe General Purpose Sling - Mesh - Extra Large 355 | 1 | | 300 kg | Mesh |