

## Stocking Aid - JOBST Stocking Donner



- A plastic coated metal stocking applicator to assist in putting on compression or surgical stockings
- The frame has foam covered handles
- The frame will hold the stocking/sock sufficiently open to allow the foot to be inserted
- The stocking still needs to be fed over the frame using the hands but may eliminate prolonged bending to reach the feet
- Available in one size and colour

### Warranty Information

- 1 Year



Code	Display Name
PTS641090	Stocking Aid - JOBST Stocking Donner

Contracted DVA product

Contract approved

**ndis**

**Item Code:**

05\_090900046\_0103\_1\_2

Click [here](#) for more information

**Compliance**

No applicable compliance standard