

Aspire Montana Single Action Maxi Lift Recline Chair



- The Montana Lift Recline Chair provides a large seating surface, offering improved comfort for users who need that little bit of extra room.
- Single motor lift recline chair
- Designed for larger users
- Head & armrest covers for improved longevity
- Superior wipe clean microsuede upholstery
- Part of the Aspire BRONZE range of lift recline chairs
- Available in Bronze Microsuede Leather-look
- Looking to protect your Electric Lift Recline Chair? Try our [Universal Chair Protector](#).

✓ Contracted DVA product

✓ Contract approved

ndis

Item Code:

05_180903180_0103_1_2

Click [here](#) for more information

Compliance

ISO13934-1 / ISO13937-1 / BS6807 / BS7175

| Code | Display Name | RRP | Box Quantity | Seat Depth | Seat Width | Backrest Height | Seat To Floor Height | Safe Working Load |
|-----------|--|------|--------------|------------|------------|-----------------|----------------------|-------------------|
| CHP227400 | Chair Powerlift Recline - Aspire Montana Bariatric Recline - Bronze - Single Motor - 220kg | 2659 | 1 | 530 mm | 630 mm | 660 mm | 510 mm | 220 Kg |