

## Forearm Crutches - Custom




- Reinforced aluminium construction designed for users with long term needs
- Absence of adjustment holes makes the crutch 'silent' in use
- Supplied in one standard length and may be cut to the required length
- Please specify handle to floor height at time of order
- Ferrule: 25 mm
- Product Weight: 800 g (each)

### Warranty Information

- 1 Year



 Contract approved

### ndis

**Item Code:**

05\_120309088\_0105\_1\_2

**Walking Aids -**

**Composite:**05\_416\_0105\_1\_2

Click [here](#) for more information

### Compliance

ISO 11199-2:2005

Code	Display Name	RRP	Box Quantity	Safe Working Load
WAC688400	Forearm Crutches - Custom	64	21	190 kg