

## MOTomed Hand Fixation



**RECK**  
MOTomed®

- Hand fixation by wrist cuff available in 3 sizes
- Easy and quick fixation of the paralysed (weak) hand on the arm/ upper body trainer, on handle or the various hand grips
- Recommended for users with limited gripping function



Code	Display Name
FPM102220	MOTomed Hand Fixation - Large - Adult
FPM102240	MOTomed Hand Fixation - Medium - Young Adult
FPM102230	MOTomed Hand Fixation - Small - Children