


Sigvaris Traditional Compression Stockings



- Sigvaris Traditional medical compression garments
- Maximum pressure stability throughout the day
- Excellent strength-stretch ratio
- Great medical benefits for safe and successful treatment
- Can be washed on a delicate wash in a washing bag and air dried
- Please specify fitting (open or closed toe), size and colour at time of ordering

- ✔ Contracted DVA product
- ✔ Contract approved

Code	Display Name	Colour
PTS650225	SIGVARIS Traditional - Knee - 23-32mmHg XL PLUS Short Beige Open Toe	Beige
PTS650155	SIGVARIS Traditional - Knee - 23-32mmHg L Long Beige Open Toe	Beige
PTS650230	SIGVARIS Traditional - Knee - 23-32mmHg XL Short Beige Open Toe	Beige
PTS650195	SIGVARIS Traditional - Knee - 23-32mmHg S Long Beige Open Toe	Beige
PTS650220	SIGVARIS Traditional - Knee - 23-32mmHg XL PLUS Long Beige Open Toe	Beige
PTS650215	SIGVARIS Traditional - Knee - 23-32mmHg XL Long Beige Open Toe	Beige
PTS650210	SIGVARIS Traditional - Knee - 23-32mmHg S PLUS Short Beige Open Toe	Beige

Code	Display Name	Colour
PTS650205	SIGVARIS Traditional - Knee - 23-32mmHg S PLUS Long Beige Open Toe	Beige
PTS650200	SIGVARIS Traditional - Knee - 23-32mmHg S Short Beige Open Toe	Beige
PTS650190	SIGVARIS Traditional - Knee - 23-32mmHg M Short Beige Open Toe	Beige
PTS650185	SIGVARIS Traditional - Knee - 23-32mmHg M PLUS Short Beige Open Toe	Beige
PTS650180	SIGVARIS Traditional - Knee - 23-32mmHg M PLUS Long Beige Open Toe	Beige
PTS650175	SIGVARIS Traditional - Knee - 23-32mmHg M Long Beige Open Toe	Beige
PTS650170	SIGVARIS Traditional - Knee - 23-32mmHg L Short Beige Open Toe	Beige
PTS650160	SIGVARIS Traditional - Knee - 23-32mmHg L PLUS Long Beige Open Toe	Beige
PTS650305	SIGVARIS Traditional - Knee w/Grip Top - 23-32mmHg S Long Beige Open Toe	Beige
PTS650235	SIGVARIS Traditional - Knee - 23-32mmHg XS Long Beige Open Toe	Beige
PTS650165	SIGVARIS Traditional - Knee - 23-32mmHg L PLUS Short Beige Open Toe	Beige
PTS650300	SIGVARIS Traditional - Knee w/Grip Top - 23-32mmHg M Norm Beige Open Toe	Beige
PTS650240	SIGVARIS Traditional - Knee - 23-32mmHg XS Short Beige Open Toe	Beige
PTS650325	SIGVARIS Traditional - Knee w/Grip Top - 23-32mmHg XS Long Beige Open Toe	Beige
PTS650295	SIGVARIS Traditional - Knee w/Grip Top - 23-32mmHg M Long Beige Open Toe	Beige
PTS650315	SIGVARIS Traditional - Knee w/Grip Top - 23-32mmHg XL Long Beige Open Toe	Beige
PTS650310	SIGVARIS Traditional - Knee w/Grip Top - 23-32mmHg S Norm Beige Open Toe	Beige

Code	Display Name	Colour
PTS650330	SIGVARIS Traditional - Knee w/Grip Top - 23-32mmHg XS Norm Beige Open Toe	Beige
PTS650320	SIGVARIS Traditional - Knee w/Grip Top - 23-32mmHg XL Norm Beige Open Toe	Beige
PTS650290	SIGVARIS Traditional - Knee w/Grip Top - 23-32mmHg L Norm Beige Open Toe	Beige
PTS650245	SIGVARIS Traditional - Knee w/Grip Top - 23-32mmHg L PLUS Long Beige Open Toe	Beige
PTS650280	SIGVARIS Traditional - Knee w/Grip Top - 23-32mmHg XL PLUS Norm Beige Open Toe	Beige
PTS650275	SIGVARIS Traditional - Knee w/Grip Top - 23-32mmHg XL PLUS Long Beige Open Toe	Beige
PTS650270	SIGVARIS Traditional - Knee w/Grip Top - 23-32mmHg S PLUS Norm Beige Open Toe	Beige
PTS650265	SIGVARIS Traditional - Knee w/Grip Top - 23-32mmHg S PLUS Long Beige Open Toe	Beige
PTS650260	SIGVARIS Traditional - Knee w/Grip Top - 23-32mmHg M PLUS Norm Beige Open Toe	Beige
PTS650285	SIGVARIS Traditional - Knee w/Grip Top - 23-32mmHg L Long Beige Open Toe	Beige
PTS650250	SIGVARIS Traditional - Knee w/Grip Top - 23-32mmHg L PLUS Norm Beige Open Toe	Beige
PTS650255	SIGVARIS Traditional - Knee w/Grip Top - 23-32mmHg M PLUS Long Beige Open Toe	Beige