

Aspire Toilet Seat Raiser with Lid





- Moulded plastic seat raiser with plastic lid to mask odours
- Fits all standard toilet seats
- Adjustable mounting brackets to allow for correct fitting to different size toilets
- Available in three different height versions to suit individual requirements
- SWL: 200 kg

Warranty Information

- 1 Year



 Contracted DVA product

 Contract approved

ndis

Item Code:

05_091218063_0103_1_2

Toilet Aids – Composite:

05_414_0103_1_2

Click [here](#) for more information

Compliance

AS/NZS 3973:2009

Code	Display Name
BTT149900	Aspire Toilet Seat Raiser with Lid- 50mm
BTT150100	Aspire Toilet Seat Raiser with Lid - 100mm
BTT150400	Aspire Toilet Seat Raiser with Lid - 150mm