

## Walking Stick Tip



- The Hugo® Claw™ Standing Cane Tip lets your cane stand on its own!
- Lets your cane stand on its own! Stability in every direction.
- 57% lighter than the QuadPod, weighing only 0.2 lb (92 grams).
- 360° grip: great for folding canes!
- Tool-free installation compatible with all canes with 3/4" (1.9 cm) cane tubing.
- There are no reviews for this product.

Contracted DVA product



**Code**

**Display Name**

WAS747429 Walking Stick Tip - Claw - Mini Stable Tip