

## Embalming / Autopsy Trolley



- 4 x buffers on corners
- Hand control
- Top section all 304 stainless steel construction
- Electric tilt function for easy draining
- Drain with hose
- Push handle both ends
- Plug in rechargeable battery pack
- 125 mm diameter total lock twin wheel castors
- Zinc treated epoxy powder coated frame
- Curved base design
- Electric high/low operation 450 - 850mm



### Warranty Information

- 3/10 Years

sova / dignity in motion

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight
FUN100300	Embalming/Autopsy Trolley - 200 kg	2200 mm	740 mm	615 - 1015 mm	200 kg	90 kg