


Medical Pulp Shapes

- A great range of medical pulp shapes to be used in macerators

 Contract approved

Compliance

ISO 15883/AS 4187:2019



Code	Display Name	Capacity
UPP100215	Slipper Pan - Midi - 1300ml	1300 ml
UPP100220	Slipper Pan - Maxi - 2000ml	2000 ml
UPP100210	Bed Pan - 2000ml	2000 ml
UPP100230	Kidney Dish - 700ml	700 ml
UPP100245	Male Urine Bottle - 900ml - Square	900 ml
UPP100250	Male Urine Bottle - 875ml - Bulbous	875 ml
UPP100235	GP Bowl - 1000ml	1000 ml

Code	Display Name	Capacity
UPP100225	Commode Bowl - 1700ml	1700 ml
UPP100240	GP Bowl - 3000ml	3000 ml