

## Solaris ExoSoft - Lymphoedema - Full Finger Glove



- This circular-knit arm sleeve provides comfortable gradient compression and moderate containment
- Wide bands at the wrist and axilla help the sleeve stay in place and resist rolling. Created to replace daytime compression bandaging, or serve as an alternative or supplement to elastic compression stockings
- Microfiber materials are comfortable and breathable while still providing containment for mild to moderate edema
- Soft and smooth materials are easier for patients to don

✓ Contracted DVA product

✓ Contract approved

### Support @ Home

**Function:** Supporting and stimulating blood circulation

**Category:** Low risk

Click [here](#) for more information

Code	Display Name	Colour
PTA537680	Solaris ExoSoft - (20-30mmHg) Full Finger Glove - Average - Beige - Medium	Beige