

## DONJOY PROCARE 79-811 SUPER LITE AFO



- Lightweight and low profile.
- Will fit inside most lace-up shoes.
- Durable, orthopaedic grade thermoplastic may be custom trimmed and reshaped. Application:
- Indicated for foot drop

- ✔ Contracted DVA product
- ✔ Contract approved

Code	Display Name
FBS305240	DONJOY PROCARE 79-811 SUPER LITE AFO - MENS - RIGHT - UNIVERSAL
FBS305250	DONJOY PROCARE 79-811 SUPER LITE AFO - MENS - LEFT - UNIVERSAL
FBS305260	DONJOY PROCARE 79-811 SUPER LITE AFO - WOMENS - RIGHT - UNIVERSAL
FBS305270	DONJOY PROCARE 79-811 SUPER LITE AFO - WOMENS - LEFT - UNIVERSAL