

## MoliCare Premium Form



- Pads for moderate to severe incontinence
- Shaped to fit perfectly to the body
- 3 layer absorbent core keeps skin dry
- Maintains the skin pH and reduces the risk of incontinence related skin irritation
- Inner leakage cuffs for maximum leakage protection
- Odour neutralisers for high level confidence
- Soft textile like materials for ultimate wearing comfort and discretion
- Recommended to be worn with MoliCare Premium Fixpants



### Compliance

No applicable compliance standard

Code	Display Name	Capacity
CON252600	MoliCare Form - Normal Plus (Pkt 30)	4 Drops (1483ml)