

Rifton Pacer Gait Trainer



The Pacer, because it can be configured in so many ways, can accommodate clients of every ability, in any environment. The Pacer is easily configured to meet your client's every need. Choose between dynamic and standard upper frames and pair them with standard, utility and treadmill/stability bases.



- Dynamic Movement: Accommodates the natural side-to-side and up and down movements of the body during step-taking to facilitate a more effective gait pattern
- Choice of Two Uppers: standard upper frame and dynamic upper frame
- Easy Height Adjustment: Both upper frames feature a single telescoping column with a gas spring to make height adjustments easy
- Versatile in Every Setting: a system built on two bases
- Quick and Easy Positioning: The Multi-position saddle (MPS) is the ultimate positioning accessory for easy transfers and gait training
- Smoother Transfers: The MPS makes transfers, especially for heavy clients, so much easier (and safer). By placing the upper frame and MPS in their lowest positions you can ease your seated clients onto the saddle.
- A Measurable Improvement: Now available on every Pacer base as an option, Rifton's new odometer makes it easy to record client progress in any setting, in either feet or meters
- Infinite Control for Effective Movement: The Pacer casters' five functions let you fine tune the Pacer's movement as your clients gain control
- Old Pacer Prompts Attach on New Frame: Prompts and supports – such as the new arm supports and old chest support – attach anywhere along the Pacer frame's top bars
- Easily Transportable: Our detachable upper frame makes it easy to stow and transport the Pacer
- Choice of Colours: The Pacer gait trainer is available in your choice of 5

colours: blue, red, lime, pink and grey

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight
PED470074	Rifton Pacer Gait Trainer (as per script)	760 - 910 mm	660 - 760 mm	550 - 690 mm	68 kg	5.2 - 6.8 kg