

# Mortuary Body Rack



- Open designs allows access for cleaning
- 3 tier modular design (multiple bays can be assembled)
- Adjustable stainless steel / Delron feet
- Customisable tier heights
- Can be built to accommodate varying size trays
- Massive 350kg per bay rated
- \*Tier Height/SWL (floor to top of rollers)
  - 1st: 477 mm – Max Weight 350 kg
  - 2nd: 1032 mm – Max Weight 300 kg
  - 3rd: 1587 mm – Max Weight 250 kg
  - 4th: 1795 mm – Max Weight 200 kg



## Warranty Information

- 10 Years

sova / dignity in motion

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight
FUN100520	3 Tier Mortuary Rack	2110 mm	740 mm	1690 mm	350 kg /Bay *	90 kg
FUN100560	4 Tier Mortuary Rack	2110 mm	740 mm	1900 mm	350 kg /Bay *	90 kg