

Freedom Viking Acute Forearm Walking Tutor - BRO214



- Designed for hospitals and other health care institutions
- Hydraulic height adjustment with easy adjustment lever
- Width adjustable soft polyurethane padded armrests
- Angle adjustable soft handgrips
- Strong aluminium tube frame supporting users up to 200 kg
- Width adjustable base
- Comes with oxygen cylinder holder

- ✔ Contracted DVA product
- ✔ Contract approved

Compliance
ISO 11199-2:2005

Warranty Information

- 1 Year



Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Material	Colour
WAF699130	Freedom Viking Acute Forearm Walking Tutor	780 mm	715 - 915 mm	920 - 1350 mm	200 kg	Aluminium	Blue
GSP650520	Spare Part - Rear Lockable Castor Only	-	-	-	-	-	-