

Oxford Advance Folding Patient Lifter



- Mobile patient lifter for transfer to a seated or lying position
- Compact design can be folded easily for storage or transportation
- Epoxy powder coated frame
- Single speed handset control
- Emergency lowering
- Over-sized push handles for excellent manoeuvrability
- Manual foot operated leg gating and tapered leg design for close access
- Suits a large range of patient slings
- Supplied with battery and charger
- Item Weight: 31.7 kg

- ✔ Contracted DVA product
- ✔ Contract approved

ndis

Item Code:

05_123603161_0105_1_2

Hoist –

Composite:05_420_0105_1_2

Click [here](#) for more information

Warranty Information

- 2 Years



Compliance

AS/NZS ISO 10535:2011

| Code | Display Name | Total Length | Total Width | Total Height | Lift Range | Leg Gait Range | Safe Working Load | Material |
|-----------|--------------------------------|--------------|-------------|----------------|---------------|----------------|-------------------|-----------|
| LSS377300 | Advance Folding Patient Lifter | 1300 mm | 550 mm | 1360 - 1860 mm | 350 - 1860 mm | 570 - 1000 mm | 155 kg | Aluminium |