

## Sphygmomanometer Digital Auto



- Fully automatic providing accurate and comfortable upper arm blood pressure measurement with Enhanced IntelliSense Technology.
- Personalised inflation for maximum comfort
- Quick deflation releases valve for speedy measurement
- Stores 30 Measurements in Memory
- Body Movement Detection
- Irregular Heartbeat Detection
- Hypertension Indicator

- ✔ Contracted DVA product
- ✔ Contract approved

### Warranty Information

- 1 Year

Code	Display Name	RRP	Box Quantity
GMD300330	Sphygmomanometer Digital Auto	234	1