

Heelift Suspension Boot



- Designed by orthopaedic surgeons specifically for prevention and treatment of heel pressure ulcers
- Effectively eliminates heel pressure by suspending the heel in the air
- Two options available: Smooth or Convolutated

✔ Contracted DVA product

✔ Contract approved

Warranty Information

- 1 Year

ndis

Item Code:

05_098800044_0103_1_2

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name	RRP	Box Quantity
PTA509200	Heelift Suspension Boot - Smooth	121	1
PTA509100	Heelift Suspension Boot - Convolutated	121	1