

Body Assist 590 Hernia Belt



- Removable foam pads
- Under leg straps
- Easy velcro fastening and adjustment
- Made of lightweight material construction & assists to inhibit skin irritation
- The belt can be worn whilst showering or swimming

- ✔ Contracted DVA product
- ✔ Contract approved

Code	Display Name	Waist/Hip Circumference
FBS3016134	BodyAssist - Hernia Aid - Hernia Belt - XL	1050 - 1150 mm
FBS301613	BodyAssist - Hernia Aid - Hernia Belt - X Sml	650 - 750 mm
FBS3016131	BodyAssist - Hernia Aid - Hernia Belt - Sml	750 - 850 mm
FBS3016132	BodyAssist - Hernia Aid - Hernia Belt - Med	850 - 950 mm
FBS3016133	BodyAssist - Hernia Aid - Hernia Belt - Lge	950 - 1050 mm