


Nosey Cup

- These cups are designed for added comfort and convenience
- A special cutout helps to maintain the proper head and neck positioning when swallowing, making drinking even easier
- Weight: 45 g

 Contract approved

ndis

Item Code:

05_150900319_0103_1_2

Click [here](#) for more information

Compliance

No applicable compliance standard



Warranty Information

- 1 Year

Code	Display Name	Capacity
DLK275480	Nosey Cup 115ml - Blue Tint	115 ml
DLK275490	Nosey Cup 225ml - Blue Tint	225 ml
DLK275500	Nosey Cup 355ml - Blue Tint	355 ml