

ROHO Cushion - High Profile - High Pressure Care



- Offers all the benefits of ROHO Dry Floatation but with less postural control and flexibility than offered by the Dual Compartment and Quadtro Select ranges
- Interconnected cells help disperse pressure and reduce shear risk promoting pressure management
- Appropriate for clients with high risk of ischemic ulcers
- Soft neoprene construction is flame retardant and durable
- Maximum User Weight: Unlimited
- High Profile (HP): 4" (10cm) High Cells, Higher pressure disbursement



*** Make and model of hire items may vary depending on availability. Minimum rental is two weeks.

Code

Display Name

H_PTC543700 ROHO Cushion - High Profile - High Pressure Care