

## Jobst Compression Stockings - Men



- Jobst men compression stockings combine style with effective compression support
- Attractive ribbed design with all-day non-constricting comfort grip band keeps socks up without binding or pinching
- Special soft micro fibres wick away moisture six times faster than standard nylon, keeping feet cooler and dryer longer
- Special formula fights odour by inhibiting growth of fungi and bacteria on the socks
- Durable, reinforced heel resists shoe abrasion, strong yarns resist wear and tear

- ✔ Contracted DVA product
- ✔ Contract approved

### Warranty Information

- 1 Year

Code	Display Name	Colour	Calf/Arm Circumference	Knee Circumference	Ankle Circumference
PTS640040	Jobst Compression Stockings - Knee - 15-20mmHg - C/Toe - Sml	Black	28 - 38 cm	40 -62 cm	18 -21 cm
PTS640530	Jobst for Men - Knee - 20-30mmHg - C/Toe - Sml	Black	28-38cm	40-62cm	18-21cm
PTS640050	Jobst Compression Stockings - Knee - 15-20mmHg - C/Toe - Med	Black	30 - 42 cm	46 -70 cm	21 - 25 cm
PTS640540	Jobst for Men - Knee - 20-30mmHg - C/Toe - Med	Black	30-42cm	46-70cm	21-25cm
PTS640060	Jobst Compression Stockings - Knee - 15-20mmHg - C/Toe - Lge	Black	32 - 46 cm	54 - 78 cm	25 - 29 cm
PTS640550	Jobst for Men - Knee - 20-30mmHg - C/Toe - Lge	Black	32-46cm	54-78cm	25-29cm

Code	Display Name	Colour	Calf/Arm Circumference	Knee Circumference	Ankle Circumference
PTS640070	Jobst Compression Stockings - Knee - 15-20mmHg - C/Toe - XL	Black	34 - 50 cm	60 - 81 cm	29 - 33 cm
PTS640560	Jobst for Men - Knee - 20-30mmHg - C/Toe - XL	Black	34-50cm	60-81cm	29-33cm
PTS640830	Jobst for Men - Knee - 30-40mmHg - C/Toe - Sml	Black	28-38cm	40-62cm	18-21cm
PTS640840	Jobst for Men - Knee - 30-40mmHg - C/Toe - Med	Black	30-42cm	46-70cm	21-25cm
PTS640850	Jobst for Men - Knee - 30-40mmHg - C/Toe - Lge	Black	32-46cm	54-78cm	25-29cm
PTS640860	Jobst for Men - Knee - 30-40mmHg - C/Toe - XL	Black	34-50cm	60-81cm	29-33cm
PTS640601	Jobst for Men - Thigh - 20-30mmHg - C/Toe - Sml	Black	28-38cm	40-62cm	18-21cm
PTS640602	Jobst for Men - Thigh - 20-30mmHg - C/Toe - Med	Black	30-42cm	46-70cm	21-25cm
PTS640603	Jobst for Men - Thigh - 20-30mmHg - C/Toe - Lge	Black	32-46cm	54-78cm	25-29cm
PTS640604	Jobst for Men - Thigh - 20-30mmHg - C/Toe - XL	Black	34-50cm	60-81cm	29-33cm
PTS640605	Jobst for Men - Thigh - 30-40mmHg - C/Toe - Sml	Black	28-38cm	40-62cm	18-21cm
PTS640606	Jobst for Men - Thigh - 30-40mmHg - C/Toe - Med	Black	30-42cm	46-70cm	21-25cm
PTS640607	Jobst for Men - Thigh - 30-40mmHg - C/Toe - Lge	Black	32-46cm	54-78cm	25-29cm
PTS640608	Jobst for Men - Thigh - 30-40mmHg - C/Toe - XL	Black	34-50cm	60-81cm	29-33cm