


## Oxford Presence Patient Lifter



- Mobile patient lifter for transfer to a seated or lying position
- High lifting range accomodates transfers to higher surfaces
- Large load capacity of 227 kg
- Robust epoxy powder coated frame construction makes it ideal for facilities
- Electric hand control operation
- Emergency lowering and raising
- Over-sized push handles for excellent manoeuvrability
- 100 mm rear locking castors and 100 mm front twin swivel castors
- Electric leg gating
- Suits a large range of patient slings
- Supplied with battery and charger
- Pictured with optional 4 point pivot frame (LSS377395)
- Item Weight: 40 kg

 Contract approved

**ndis**

**Item Code:**

05\_123603161\_0105\_1\_2

**Hoist –**

**Composite:**05\_420\_0105\_1\_2

Click [here](#) for more information

**Compliance**

AS/NZS ISO 10535:2011



### Warranty Information

- 2 Years

Code	Display Name	Total Length	Total Width	Total Height	Lift Range	Leg Gait Range	Electric Leg	Safe Working Load	Material
LSS377325	Presence Patient Lifter	1460 mm	670 mm	1460 - 2061 mm	440 - 1910 mm	550 - 1020 mm	True	227 kg	Aluminium