


Pride 660 Powerlift Recliner - Dual Motor



- To suit the bariatric user
- Engineered furniture laminate/hardwood frame
- Patented, quiet & smooth lift system
- Large side pocket
- Hand control can be easily moved to left-hand side
- Integrated emergency battery backup (batteries not included)
- External transformer for easy access to battery backup

 Contract approved

ndis

Item Code:

05_180903180_0103_1_2

Click [here](#) for more information

Warranty Information

- 2 Years

Compliance

ISO13934-1 / ISO13937-1 / BS6807
/ BS7175

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

Code	Display Name	Seat Depth	Seat Width	Backrest Height	Seat To Floor Height	User Height	Safe Working Load	Colour
CHP193702	660 Recliner - Latte	510 mm	510 mm	635 mm	480 mm	< 1630 mm	170 kg	Latte Fabric