

Liko Amputee Sling (Mod. 75)




- Designed for safe and secure lifting of patients with amputated legs
- Supports patient with single and double leg amputations as well as patients with contractures
- Can also be used for non-amputee patients who are slender and require a smaller sling
- Offers a sitting posture that is slightly reclined with leg supports applied parallel to each under both thighs
- Model 75 offers head support and is available in Polyester

Warranty Information

- 1 Year



 Contract approved

ndis

Item Code:

05_123603166_0105_1_2

Hoist –

Composite:05_420_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Safe Working Load	Material
LSS394520	Amputee Highback Sling - Reinforced Leg Support - Polyester - SMALL	200 kg	Polyester (Loops)
LSS394530	Amputee Highback Sling - Reinforced Leg Support - Polyester - MEDIUM	200 kg	Polyester (Loops)
LSS394540	Amputee Highback Sling - Reinforced leg Support - Polyester - LRG	200 kg	Polyester (Loops)
LSS394550	Amputee Highback Sling - Reinforced Leg Support - Polyester - XLARGE	300 kg	Polyester (Loops)