

Rehband 1068 EQ Posture Support



- A flexible aid for training a better posture
- Bendable stays for additional support
- Elastic and well-fitting
- Wear with singlet or t-shirt underneath
- Sizing: waist circ (cm)

✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Upper limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Chest Circumference
FBS301258	Rehband - Shoulder Support - EQ posture support - X Sml	560 - 640 mm
FBS301258L	Rehband - Shoulder Support - EQ posture support - Lge	800 - 900 mm
FBS301258M	Rehband - Shoulder Support - EQ posture support - Med	720 - 800 mm
FBS301258S	Rehband - Shoulder Support - EQ posture support - Sml	640 - 720 mm
FBS301258XL	Rehband - Shoulder Support - EQ posture support - XL	900 -1000 mm