

Inner Lip Plate



- Plate with an inner lip to hold food on plate while user scoops against it
- Polycarbonate sandstone colour
- Not for use with serrated knives
- Dishwasher and microwave safe up to 80°C
- Weight: 136 g

✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Eating and drinking

Category: Low risk

Click [here](#) for more information

NDIS

Item Code:

05_150900319_0103_1_2

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name	External Diameter
------	--------------	-------------------

DLK275290 Inner Lip Plate 228 mm