

## Aidacare Medication Cart



- The Aidacare Medication Cart is the powerful all-round solution for a multitude of applications including Treatment, Procedure, IV, Trauma, Dressing, Supply and General cart
- Aidacare carts are designed to offer an attractive, practical and durable cart with multiple configurations to enable the efficient organisation and delivery of medication and equipment
- Robust panels to minimise dents and scratches
- Large work area with slide out work surface are standard
- Easy roll 125 mm castors
- Fully customisable with a host of accessories and drawer configurations.
- **Standard Inclusions:** Drawer Dividers and Extended Sliding Surface
- Lock System Options:
  - **Key lock:** The key lock system is a proven method of securing your medical cart and providing the level of security necessary
  - **Digital lock:** A user friendly system where the cart is locked and unlocked using a pre-set code ensuring ease of access and the security required in clinical environments. The digital lock is robust and waterproof
- Standard drawer configuration detailed below
- **Small Medication Cart:**
  - 1 small 3' drawer (7 cm)
  - 2 medium 6' drawers (16 cm)
  - 1 large 9' drawer (25 cm)
- **Medium Medication Cart:**
  - 2 small 3' drawers (7 cm)
  - 2 medium 6' drawers (16 cm)
  - 1 large 9' drawer (25 cm)
- **Large Medication Cart:**



### Accessories

Click on the cart icon to display the range of optional accessories

- 1 small 3' drawer (7 cm)
- 3 medium 6' drawers (16 cm)
- 1 large 9' drawer (25 cm)

## Warranty Information

- Structural: 5 Years Warranty
- Digital Lock: 2 Years Warranty
- Accessories: 1 Year Warranty

| Code      | Display Name   | Total Length | Total Width | Total Height | Colour |
|-----------|--|--------------|-------------|--------------|--------|
| WAT909430 | Aidacare Medication Cart - Small - Digital Lock - Blue   | 790 mm       | 610 mm      | 900 mm       | Blue   |
| WAT909510 | Aidacare Medication Cart - Small - Digital Lock - Cream  | 790 mm       | 610 mm      | 900 mm       | Cream  |
| WAT909470 | Aidacare Medication Cart - Small - Key Lock - Blue       | 790 mm       | 610 mm      | 900 mm       | Blue   |
| WAT909560 | Aidacare Medication Cart - Small - Key Lock - Cream      | 790 mm       | 610 mm      | 900 mm       | Cream  |
| WAT909450 | Aidacare Medication Cart - Medium - Digital Lock - Blue  | 790 mm       | 610 mm      | 1000 mm      | Blue   |
| WAT909540 | Aidacare Medication Cart - Medium - Digital Lock - Cream | 790 mm       | 610 mm      | 1000 mm      | Cream  |
| WAT909490 | Aidacare Medication Cart - Medium - Key Lock - Blue      | 790 mm       | 610 mm      | 1000 mm      | Blue   |
| WAT909590 | Aidacare Medication Cart - Medium - Key Lock - Cream     | 790 mm       | 610 mm      | 1000 mm      | Cream  |
| WAT909460 | Aidacare Medication Cart - Large - Digital Lock - Blue   | 790 mm       | 610 mm      | 1100 mm      | Blue   |
| WAT909550 | Aidacare Medication Cart - Large - Digital Lock - Cream  | 790 mm       | 610 mm      | 1100 mm      | Cream  |

| Code      | Display Name  | Total Length | Total Width | Total Height | Colour |
|-----------|---|--------------|-------------|--------------|--------|
| WAT909500 | Aidacare Medication Cart - Large - Key Lock - Blue  | 790 mm       | 610 mm      | 1100 mm      | Blue   |
| WAT909600 | Aidacare Medication Cart - Large - Key Lock - Cream | 790 mm       | 610 mm      | 1100 mm      | Cream  |