

## JAY - Care Backrest



**JAY**<sup>®</sup>

- Designed to accommodate a client who has a fixed kyphotic, the JAY-CARE backrest provides solid back support to accommodate and slow-down back deformity
- The contoured design and 39° angle adjustability help further accommodate the multiple needs of the kyphotic client and make the JAY-CARE suitable to a client who has changing postural needs over time
- Provides mild to moderate lateral support through the built-in lateral contour of the shell
- Optional foam wedges provide for a more customised fit
- All models include: Jay AirExchange cover and universal mounting hardware

### Warranty Information

- 2 Years Backrest
- 6 Months Cover



**ndis**

**Item Code:**

05\_181003192\_0105\_1\_2

**Wheelchair Manual – Composite:**

**Wheelchair Powered –**

**Composite:**05\_417\_0105\_1\_2

05\_419\_0105\_1\_2

Click [here](#) for more information

### Compliance

No applicable compliance standard

Code	Display Name	Seat Width	Backrest Width	Backrest Height	Safe Working Load
PTB532530	JAY Care - 16in	15 - 16 in	16 in	16.5 in	113 kg
PTB532520	JAY Care - 18in	17 - 18 in	18 in	17.25 in	113 kg
PTB532515	JAY Care - 20in	19 - 20 in	20 in	19 in	113 kg