

Elevating Legrest



- A great range of accessories to fit the Aspire Manual Wheelchair range to improve functionality and user experience
- Click [here](#) to ensure this accessory fits your specific wheelchair model

- ✔ Contracted DVA product
- ✔ Contract approved

Warranty Information

- 1 Year

Compliance

AS3696-1:2008



Code	Display Name	RRP	Box Quantity
WCA848100	Elevating Legrest - Left - 300/320/350/400/410mm	198	1
WCA848110	Elevating Legrest - Right - 300/320/350/400/410mm	198	1
WCA848120	Elevating Legrest - Left - 450/460/500/510mm	198	1
WCA848130	Elevating Legrest - Right - 450/460/500/510mm	198	1
WCA848140	Elevating Legrest - Left - 570/610mm	198	1
WCA848150	Elevating Legrest - Right - 570/610mm	198	1