

## MOTomed Loop.LA



**RECK**  
MOTomed®

The MOTomed Loop brings movement therapy into the future. Through the new design and its digital intelligence the user can experience passive or active training in a completely new way. With the new and innovative height adjustment, users can adjust the position of the MOTomed loop without using tools. The Loop LA features full software functionality, including disability-specific training programs and some interactive games



- The LA comes with both upper limb and lower limb trainers installed
- PVC coated safety foot shells (pair, incl. pedal cranks)
- Soft Grip Handles (pair, incl. pedal cranks)
- Gas spring support for height adjustment
- Special stand with reduced stand height
- Leg guides with plastic-coated calf shells (pair)
- QuickFix foot fastening system (pair)
- Hand fixation by wrist cuff (piece)
- Serial port activated: trainer can be upgraded to an FES Cycle with the HASOMED RehaMove2 'Functional Electrical Stimulator'

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight
FPM101130	MOTomed Loop.LA - Serial Interface, Gas Spring, Leg guides, Foot & Hand Fix	70 cm	60 cm	107 - 122 cm	135 kg	36 kg