

## Raised Toilet Seat with Swing Back Arms



- This raised toilet seat features swing back arms and is height adjustable
- Fits most standard toilets
- Fixes onto the toilet bowl for safe transferring (replaces existing seat & lid.
- Plastic moulded armrests
- Height adjustable range: 5, 10 and 15 cm
- Width Between Arms: 410 mm

### Warranty Information

- 1 Year

Contracted DVA product

Contract approved

### ndis

#### Item Code:

05\_091218063\_0103\_1\_2

#### Toilet Aids – Composite:

05\_414\_0103\_1\_2

Click [here](#) for more information

### Compliance

AS/NZS 3973:2009

Code	Display Name
BTT150450	Raised Toilet Seat with Swing Back Arms