

Economy Clean Linen Trolley - Open



- Open both sides for easy access
- Mobile clean linen trolley
- Transport clean linen in the hospital or aged care environment
- Lightweight construction
- 125 mm castors, wall buffer discs, two shelves, top tray and push handle



Compliance

No applicable compliance standard

Warranty Information

- 3 Years

Code	Display Name	Total Length	Total Width	Total Height
LAC353640	Clean Linen Trolley - Small	920 mm	525 mm	1040 mm
LAC353610	Clean Linen Trolley - Large	1100 mm	550 mm	1400 mm