

## K Care Commode Seat - Wedge



- Provides specific pelvic positioning
- Marine Grade vinyl is wrapped, stapled and sealed, making cleaning easier with no seam
- Open and closed configuration available
- Comes with plastic seat runners for easy attachment and removal for cleaning



**ndis**

**Item Code:**

05\_093330084\_0103\_1\_2

**Mobile Shower Commode –**

**Composite:** 05\_413\_0103\_1\_2

Click [here](#) for more information

### Warranty Information

- 1 Year

**Compliance**

AS/NZS 3973:2009

Code	Display Name	Seat Depth	Seat Width
BTC074100	Closed Front Commode Seat - Hand Made - Wedge	430 mm	440 mm