


## Liko Viking L Lifter



- An everyday flexible lifting solution
- Informative control box displays number of lift cycles, battery level, lift overload and service indicator
- Optional quick-release hook system which means caregivers can easily and rapidly change a slingbar or install a scale
- Optional armrests help support patients during ambulation as well as making it easier for caregivers to manoeuvre the lift
- Three variable height settings accommodate a variety of different lifting situations
- Intuitive and ergonomically-designed hand control incorporates indicators for low battery level and when the Safe Working Load of the lift is exceeded
- An extensive assortment of accessories is available including slings, slingbars and scales
- Lifting Height 460 - 1790 mm

 Contract approved

**ndis**

**Item Code:**

05\_123603161\_0105\_1\_2

**Hoist –**

**Composite:**05\_420\_0105\_1\_2

Click [here](#) for more information

**Compliance**

AS/NZS ISO 10535:2011

### Warranty Information

- Frames: 10 Years
- Electrical Parts: 2 Years

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Material
LSS393180	Liko Viking Lifter - LARGE	1360 mm	695 - 1070 mm	1390 - 2030 mm	250 kg	Aluminium
LSS393210	Liko Viking L - Leg Protector Set -	-	-	-	-	-