

Jobst Relief Compression Stockings - Unisex



- Jobst Relief stockings with variety of sizes, lengths and compression levels
- Relief is an economical, unisex compression stocking
- 3-D knit structure using air covered spandex yarn that is soft to touch and comfortable for all day wear
- Reinforced heel, designed for long lasting durability
- Relief Chaps available in left and right leg open toe only with wide comfortable adjustable belt

- ✔ Contracted DVA product
- ✔ Contract approved

Warranty Information

- 1 Year

Code	Display Name	Colour	Calf/Arm Circumference	Waist/Hip Circumference	Knee Circumference	Ankle Circumference
PTS642516	Jobst Relief Compression Stockings - Knee - 15-20mmHg - O/Toe - Sml	Beige	28 - 38 cm	-	-	18 -21 cm
PTS642520	Jobst Relief - Knee - 15-20mmHg - O/Toe - Lge-Full Calf - Beige	Beige	28-38cm	-	-	18-21cm
PTS642517	Jobst Relief Compression Stockings - Knee - 15-20mmHg - O/Toe - Medium	Beige	30 - 42 cm	-	-	21 - 25 cm
PTS642521	Jobst Relief - Knee - 15-20mmHg - O/Toe - XL/Full Calf -	Beige	30-42cm	-	-	21-25cm

Code	Display Name	Colour	Calf/Arm Circumference	Waist/Hip Circumference	Knee Circumference	Ankle Circumference
	Beige					
PTS642518	Jobst Relief Compression Stockings - Knee - 15-20mmHg - O/Toe - Lge	Beige	32 - 46 cm	-	-	25 - 29 cm
PTS642526	Jobst Relief - Knee - 15-20mmHg - C/Toe - Lge/Full Calf - Beige	Beige	32-46cm	-	-	25-29cm
PTS642519	Jobst Relief Compression Stockings - Knee - 15-20mmHg - O/Toe - XL	Beige	34 - 50 cm	-	-	29 - 33 cm
PTS642527	Jobst Relief - Knee - 15-20mmHg - C/Toe - XL/Full Calf - Beige	Beige	34-50cm	-	-	29-33cm
PTS641110	Jobst Relief - Knee - 20-30mmHg - C/Toe - Sml	Beige	28-38cm	-	-	18-21cm
PTS641120	Jobst Relief - Knee - 20-30mmHg - C/Toe - Med	Beige	30-42cm	-	-	21-25cm
PTS642522	Jobst Relief Compression Stockings - Knee - 15-20mmHg - C/Toe - Sml	Beige	28 - 38 cm	-	-	18 -21 cm
PTS641130	Jobst Relief - Knee - 20-30mmHg - C/Toe - Lge	Beige	32-46cm	-	-	25-29cm
PTS641140	Jobst Relief - Knee - 20-30mmHg - C/Toe - XL	Beige	34-50cm	-	-	29-33cm
PTS642523	Jobst Relief Compression Stockings - Knee - 15-20mmHg - C/Toe - Medium	Beige	30 - 42 cm	-	-	21 - 25 cm
PTS641141	Jobst Relief - Knee - 20-30mmHg - C/Toe - Lge/Full Calf	Beige	28-38cm	-	-	18-21cm
PTS641142	Jobst Relief - Knee - 20-30mmHg - C/Toe - XL/Full Calf	Beige	30-42cm	-	-	21-25cm
PTS642524	Jobst Relief Compression Stockings - Knee - 15-20mmHg - C/	Beige	32 - 46 cm	-	-	25 - 29 cm

Code	Display Name	Colour	Calf/Arm Circumference	Waist/Hip Circumference	Knee Circumference	Ankle Circumference
	Toe - Lge					
PTS641150	Jobst Relief - Knee - 20-30mmHg - O/Toe - Sml	Beige	32-46cm	-	-	25-29cm
PTS641160	Jobst Relief - Knee - 20-30mmHg - O/Toe - Med	Beige	34-50cm	-	-	29-33cm
PTS642525	Jobst Relief Compression Stockings - Knee - 15-20mmHg - C/Toe - XL	Beige	34 - 50 cm	-	-	29 - 33 cm
PTS641170	Jobst Relief - Knee - 20-30mmHg - O/Toe - Lge	Beige	28-38cm	-	-	18-21cm
PTS641181	Jobst Relief - Knee - 20-30mmHg - O/Toe - Lge/Full Calf	Beige	32-46cm	-	-	25-29cm
PTS641180	Jobst Relief - Knee - 20-30mmHg - O/Toe - XL	Beige	30-42cm	-	-	21-25cm
PTS641182	Jobst Relief - Knee - 20-30mmHg - O/Toe - XL/Full Calf	Beige	34-50cm	-	-	29-33cm
PTS641190	Jobst Relief - Knee - 30-40mmHg - O/Toe - Sml	Beige	28-38cm	-	-	18-21cm
PTS641200	Jobst Relief - Knee - 30-40mmHg - C/Toe - Med	Beige	30-42cm	-	-	21-25cm
PTS641210	Jobst Relief - Knee - 30-40mmHg - C/Toe - Lge	Beige	32-46cm	-	-	25-29cm
PTS641221	Jobst Relief - Knee - 30-40mmHg - O/Toe - Lge/Full Calf	Beige	28-38cm	-	-	18-21cm
PTS641220	Jobst Relief - Knee - 30-40mmHg - C/Toe - XL	Beige	34-50cm	-	-	29-33cm
PTS641222	Jobst Relief - Knee - 30-40mmHg - O/Toe - XL/Full Calf	Beige	30-42cm	-	-	21-25cm
PTS641230	Jobst Relief - Knee - 30-40mmHg - O/Toe - Sml	Beige	32-46cm	-	-	25-29cm
PTS641240	Jobst Relief - Knee - 30-40mmHg - O/Toe - Med	Beige	34-50cm	-	-	29-33cm
PTS641250	Jobst Relief - Knee - 30-40mmHg - O/Toe - Lge	Beige	28-38cm	-	-	18-21cm

Code	Display Name	Colour	Calf/Arm Circumference	Waist/Hip Circumference	Knee Circumference	Ankle Circumference
PTS641261	Jobst Relief - Knee - 30-40mmHg - C/Toe - Lge/Full Calf	Beige	32-46cm	-	-	25-29cm
PTS641260	Jobst Relief - Knee - 30-40mmHg - O/Toe - XL	Beige	30-42cm	-	-	21-25cm
PTS641262	Jobst Relief - Knee - 30-40mmHg - C/Toe - XL/Full Calf	Beige	34-50cm	-	-	29-33cm
PTS641265	Jobst Relief - Thigh - 20-30mmHg - C/Toe - Sml	Beige	28-38cm	-	40-62cm	18-21cm
PTS641270	Jobst Relief - Thigh - 20-30mmHg - C/Toe - Med	Beige	30-42cm	-	46-70cm	21-25cm
PTS641280	Jobst Relief - Thigh - 20-30mmHg - C/Toe - Lge	Beige	32-46cm	-	54-78cm	25-29cm
PTS641290	Jobst Relief - Thigh - 20-30mmHg - O/Toe - XL	Beige	34-50cm	-	60-81cm	29-33cm
PTS641300	Jobst Relief - Thigh - 20-30mmHg - O/Toe - Sml	Beige	28-38cm	-	40-62cm	18-21cm
PTS641310	Jobst Relief - Thigh - 20-30mmHg - O/Toe - Med	Beige	30-42cm	-	46-70cm	21-25cm
PTS641319	Jobst Relief - Thigh - 20-30mmHg - O/Toe - Lge	Beige	32-46cm	-	54-78cm	25-29cm
PTS641320	Jobst Relief - Thigh - 20-30mmHg - C/Toe - XL	Beige	34-50cm	-	60-81cm	29-33cm
PTS641940	Jobst Relief - Thigh - 30-40mmHg - O/Toe - Sml	Beige	28-38cm	-	40-62cm	18-21cm
PTS641941	Jobst Relief - Thigh - 30-40mmHg - O/Toe - Med	Beige	30-42cm	-	46-70cm	21-25cm
PTS641942	Jobst Relief - Thigh - 30-40mmHg - O/Toe - Lge	Beige	32-46cm	-	54-78cm	25-29cm
PTS641943	Jobst Relief - Thigh - 30-40mmHg - O/Toe - XL	Beige	34-50cm	-	60-81cm	29-33cm
PTS641944	Jobst Relief - Thigh - 30-40mmHg - C/Toe - Sml	Beige	28-38cm	-	40-62cm	18-21cm
PTS641945	Jobst Relief - Thigh - 30-40mmHg - C/Toe - Med	Beige	30-42cm	-	46-70cm	21-25cm
PTS641946	Jobst Relief - Thigh - 30-40mmHg - C/Toe - Lge	Beige	32-46cm	-	54-78cm	25-29cm

Code	Display Name	Colour	Calf/Arm Circumference	Waist/Hip Circumference	Knee Circumference	Ankle Circumference
PTS642264	Jobst Relief - Thigh - 30-40mmHg - C/Toe - XL	Beige	34-50cm	-	60-81cm	29-33cm
PTS641938	Jobst Relief - Waist - 20-30mmHg - C/Toe - Sml	Beige	28-38cm	71-117cm	40-62cm	18-21cm
PTS641330	Jobst Relief - Waist - 20-30mmHg - C/Toe - Med	Beige	30-42cm	76-127cm	46-70cm	21-25cm
PTS641939	Jobst Relief - Waist - 20-30mmHg - C/Toe - Lge	Beige	32-46cm	81-137cm	54-78cm	25-29cm
PTS641335	Jobst Relief - Waist - 20-30mmHg - C/Toe - XL	Beige	34-50cm	102-166cm	60-81cm	29-33cm
PTS641339	Jobst Relief - Waist - 20-30mmHg - O/Toe - Sml	Beige	28-38cm	71-117cm	40-62cm	18-21cm
PTS641340	Jobst Relief - Waist - 20-30mmHg - O/Toe - Med	Beige	30-42cm	76-127cm	46-70cm	21-25cm
PTS641341	Jobst Relief - Waist - 20-30mmHg - O/Toe - Lge	Beige	32-46cm	81-137cm	54-78cm	25-29cm
PTS641342	Jobst Relief - Waist - 20-30mmHg - O/Toe - XL	Beige	34-50cm	102-166cm	60-81cm	29-33cm
PTS641947	Jobst Relief - Waist - 30-40mmHg - O/Toe - Sml	Beige	28-38cm	71-117cm	40-62cm	18-21cm
PTS641948	Jobst Relief - Waist - 30-40mmHg - O/Toe - Med	Beige	30-42cm	76-127cm	46-70cm	21-25cm
PTS641949	Jobst Relief - Waist - 30-40mmHg - O/Toe - Lge	Beige	32-46cm	81-137cm	54-78cm	25-29cm
PTS641950	Jobst Relief - Waist - 30-40mmHg - O/Toe - XL	Beige	34-50cm	102-166cm	60-81cm	29-33cm
PTS641349	Jobst Relief - Chaps - Left - 20-30mmHg - Sml	Beige	28-38cm	71-117cm	40-62cm	18-21cm
PTS641351	Jobst Relief - Chaps - Left - 20-30mmHg - Med	Beige	30-42cm	76-127cm	46-70cm	21-25cm
PTS641352	Jobst Relief - Chaps - Left - 20-30mmHg - Lge	Beige	32-46cm	81-137cm	54-78cm	25-29cm
PTS641350	Jobst Relief - Chaps - Left - 20-30mmHg - XL	Beige	34-50cm	102-166cm	60-81cm	29-33cm

Code	Display Name	Colour	Calf/Arm Circumference	Waist/Hip Circumference	Knee Circumference	Ankle Circumference
PTS641359	Jobst Relief - Chaps - Right - 20-30mmHg - Sml	Beige	28-38cm	71-117cm	40-62cm	18-21cm
PTS641360	Jobst Relief - Chaps - Right - 20-30mmHg - Med	Beige	30-42cm	76-127cm	46-70cm	21-25cm
PTS641365	Jobst Relief - Chaps - Right - 20-30mmHg - Lge	Beige	32-46cm	81-137cm	54-78cm	25-29cm
PTS641370	Jobst Relief - Chaps - Right - 20-30mmHg - XL	Beige	34-50cm	102-166cm	60-81cm	29-33cm
PTS641371	Jobst Relief - Chaps - Left - 30-40mmHg - Sml	Beige	28-38cm	71-117cm	40-62cm	18-21cm
PTS641380	Jobst Relief - Chaps - Left - 30-40mmHg - Med	Beige	30-42cm	76-127cm	46-70cm	21-25cm
PTS641372	Jobst Relief - Chaps - Left - 30-40mmHg - Lge	Beige	32-46cm	81-137cm	54-78cm	25-29cm
PTS641373	Jobst Relief - Chaps - Left - 30-40mmHg - XL	Beige	34-50cm	102-166cm	60-81cm	29-33cm
PTS641374	Jobst Relief - Chaps - Right - 30-40mmHg - Sml	Beige	28-38cm	71-117cm	40-62cm	18-21cm
PTS641375	Jobst Relief - Chaps - Right - 30-40mmHg - Med	Beige	30-42cm	76-127cm	46-70cm	21-25cm
PTS641376	Jobst Relief - Chaps - Right - 30-40mmHg - Lge	Beige	32-46cm	81-137cm	54-78cm	25-29cm
PTS641377	Jobst Relief - Chaps - Right - 30-40mmHg - XL	Beige	34-50cm	102-166cm	60-81cm	29-33cm