

## Bariatric Bed



- Electronic High/Low with Backrest, Knee Break and Trendelenburg
- Safe low voltage power supply
- Anti-jam backrest & Knee Break motor
- 10 button soft touch pendent handset with auto profiling
- Removable moisture resistant MDF head & foot panels
- Perforated sheet metal mattress deck for easy cleaning
- Enhanced infection control & mattress breathability
- Epoxy coated & anti rust treated steel frame
- Teflon coated pivot bushes for smooth & quiet operation
- 4 x 100 mm central & directional lockable castors
- 2 overhead pole sockets
- Under bed clearance: 150 mm



### ndis

**Item Code:**

05\_181210200\_0103\_1\_2

**Power Bed – Composite:**

05\_423\_0103\_1\_2

Click [here](#) for more information

### Compliance

IEC 60601-2-52:2009

## Warranty Information

- 15 Years Frame
- 5 Years Mechanical
- 3 Years Electrics

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Number of Functions
BEB045085	Bariatric Bed - 4 Function - King Single - 400kg	2120 mm	1050 mm	260 - 750 mm	400 kg	4