


## Walking Tutor - Adult



- Provides assistance for those who have difficulty walking.
- Tutors the user into an upright gait.
- Folding design for storage or transportation.
- Height adjustable forearm supports to suit user requirements.
- Adjustable hand grips to suit user forearm length.
- Four swivel castors for manoeuvrability.
- Powder coated mild steel construction.
- Skis or glides available as an option as well as padded arm rest.
- Small size designed for children.

### Warranty Information

- 1 Year

 Contract approved

**ndis**

**Item Code:**

05\_120606094\_0105\_1\_2

**Walking Aids -**

**Composite:**05\_416\_0105\_1\_2

Click [here](#) for more information

**Compliance**

ISO 11199-2:2005

Code	Display Name	Total Height	Safe Working Load
WAF708900	Walking Tutor - Adult	850 - 1350 mm	100 kg
WAF709000	Walking Tutor - Adult w/ Glides	850 - 1350 mm	100 kg
WAF709100	Walking Tutor - Adult w/ Padded Chest Rest	850 - 1350 mm	100 kg