

Caring Cutlery



- A range of stainless steel utensils moulded into contoured, ivory-coloured plastic handles with built-up ends. These are of value to those with a weak grip and restricted wrist or finger movement
- The knife and fork have a shaped indent on the top for the index finger, to help with directional control
- The knife also has a serrated blade to assist with cutting. In addition to the standard spoon, there is a right handed angled spoon with an angled handle to minimise wrist movement
- Hygienically sealed, easy to clean and dishwasher safe
- Handle Length: 127 mm
- Weight: 50 - 65 g

Warranty Information

- 1 Year

✔ Contracted DVA product

✔ Contract approved

ndis

Item Code:

05_150900319_0103_1_2

Click [here](#) for more information

Compliance

No applicable compliance standard

| Code | Display Name |
|-----------|--------------------------------|
| DLK276300 | Knife |
| DLK276310 | Fork |
| DLK276320 | Spoon |
| DLK276330 | Caring Angled Spoon Right Hand |