

## Chelsea Bedside Cabinet



- Simple, understated styling is captured in the Chelsea Bedside Cabinet.
- Bottom compartment features a height adjustable internal shelf for greater storage flexibility.
- Lockable top drawer keeps essential items safe and contained.
- Raised safety edging on the cabinet top prevents items from sliding off.
- Slimline handles and smooth drawer action offer function and practicality.
- Rear towel rail offers a practical and tidy solution for compact living.
- A versatile bedside cabinet to suit most settings and decor.
- Optional hutch provides added shelf space along with practical storage for books and magazines.
- Laminate finish is durable and easy to care for.
- Select from our range of laminate colour options.



Code	Display Name	Total Length	Total Width	Total Height
BEA826320	Chelsea Bedside Cabinet - Laminate	450 mm	465 mm	810 mm
BEA826330	Chelsea Bedside Cabinet w/Hutch - Laminate	450 mm	465 mm	1020 mm