


## Aspire High Back Classic Day Chair



- A high back orthopaedic chair that offers adjustability, comfort and support
- Contoured Backrest – High back design provides increased postural support
- Softly padded armrests, seat and back upholstery are comfortable and supportive
- Easy clean vinyl surfaces
- Height Adjustable – Easy push button adjustment for comfortable positioning and changing needs
- Sturdy Design - Generous 160 kg weight capacity, yet lightweight for easy manoeuvrability
- Powder Coated Steel Frame – Easy to clean, strong and durable
- Non-Slip Tips – Provide stability whilst sitting or transferring
- Available in either Slate or Champagne Vinyl

 Contracted DVA product

 Contract approved

### ndis

#### Item Code:

05\_188800051\_0103\_1\_2

Click [here](#) for more information

### Compliance

AS/NZS 4688.2:2000; AS/NZS 4688.3:2000

## Warranty Information

- 1 Year

Code	Display Name	RRP	Box Quantity	Seat Depth	Seat Width	Backrest Height	Seat To Floor Height	Safe Working Load	Colour
CHP208050	Aspire High Back Classic Day Chair - Champagne Vinyl	482	1	490 mm	490 mm	720 mm	425 - 610 mm	160 kg	Champagne Vinyl
CHP208055	Aspire High Back Classic Day Chair - Slate Vinyl	482	1	490 mm	490 mm	720 mm	425 - 610 mm	160 kg	Slate Vinyl