

## Air-Lo 8 Alternating Mattress Replacement - 200mm



The Air-Lo 8 Alternating Mattress Replacement leads the field in affordable prevention of pressure areas and skin breakdowns. It's compact, powerful yet almost silent control unit can quickly support patient loads of up to 200 kg whilst alternately allowing pressure to be released across the rest of the body.



- Static head section eliminates patient dizziness and disorientation
- 20 pure TPU transverse air cells alternate to provide maximum comfort prevention and therapy
- High quality two way stretch top cover is fully welded and can be high temperature washed for total infection control
- All waterfalls are turned and welded for maximum infection control and zip protection
- High capacity yet unobtrusive side CPR release gives maximum safety in the event of pulmonary emergency
- Medical Grade Power cord is run inside the mattress to reduce trip hazards

\*\*\* Make and model of hire items may vary depending on availability. Minimum rental is two weeks.

**Code**

**Display Name**

H\_PT636635 Air-Lo 8 Alternating Mattress Replacement - Single