

Freedom Forearm Walker - BRO212B-HD



- Swivel adjustable handles
- Padded vinyl cover for user safety
- Press button directional height adjustment
- Rear frame widened for maximum walking space
- Large indoor/ outdoor castors fitted in a U-shaped fork
- 4 Swivel Castors
- Heavy duty aluminium oval tube frame supports users up to 220 kg
- Replacement padded armrest (WAA687300)



Compliance

ISO 11199-2:2005

Warranty Information

- 1 Year



Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Material	Unit Weight	Colour
WAF699201R	Freedom Forearm Walker - Red	880 mm	660 mm	1030mm - 1240mm	220 kg	Aluminium	9.3 kg	Red